

# ULTRASONIC INTERFEROMETER (FOR LIQUIDS)

Model BL02



**Simple Instrument to measure velocity of sound in liquids**

● Working Principle

The Principle used in measuring the velocity by means of an Ultrasonic Interferometer is simple and is based on the accurate determination of wavelength ( $\lambda$ ) in the medium. If the frequency of sound vibration ( $f$ ) is accurately known, the sound velocity ( $V$ ) can be calculated according to the relation

Velocity = wavelength  $\times$  frequency

$$V = \lambda \times f$$



**Apparatus consists of:**

- High Frequency Generator
- Measuring Cell
- Connecting Cable

**Salient Features:**

- Very Small amount of liquid required
- Accuracy of measurement depends upon Micrometer provided
- Available with Micrometer of 0.01mm and 0.001mm Least Count
- Can be used for measuring over a wide range of temperature
- Jacket for water circulation to maintain constant temperature
- Measurement can be made on liquids
- No danger of change like depolymerisation in liquid as very low ultrasonic energy is used for measurement

**Improved Circuit -  
Easy Adjustment of  
Gain & Meter**

